

UGT2B4

Reactivity:Human

Tested applications:WB

Recommended Dilution:WB 1:200 - 1:2000

Calculated MW:61kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human UGT2B4

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

HLUG25; UDPGTH1; UGT2B11; UDPGT2B4; UDPGTh-1;

To place an order, please [Click HERE](#).

Catalog #:A4180

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:7363

Isotype:IgG

Swiss Prot:P06133

Purity:Affinity purification

For research use only.