

REEP5

Reactivity:Human Mouse

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:21kDa

Observed MW:Refer to figures

Immunogen:

A synthetic peptide of human REEP5

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

DP1; TB2; YOP1; D5S346; C5orf18;

To place an order, please [Click HERE](#).

Catalog #:A4191

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:7905

Isotype:IgG

Swiss Prot:Q00765

Purity:Affinity purification

For research use only.