

UXT

Reactivity:Human

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:20kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human UXT

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

STAP1; ART-27;

Catalog #:A4217

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:8409

Isotype:IgG

Swiss Prot:Q9UBK9

Purity:Affinity purification

For research use only.

Background:

The protein encoded by this gene functions as a cofactor that modulates androgen receptor-dependent transcription, and also plays a critical role in tumor necrosis factor-induced apoptosis. Expression of this gene may play a role in tumorigenesis. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene.

To place an order, please [Click HERE](#).