

## MAP4K3

**Reactivity:** Human Mouse Rat

**Tested applications:** WB IHC

**Recommended Dilution:** WB 1:500 - 1:2000 IHC 1:20 - 1:200

**Calculated MW:** 101kDa

**Observed MW:** Refer to figures

**Immunogen:**

A synthetic Peptide of human MAP4K3

**Storage Buffer:**

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

GLK; MEKKK3; MEKKK 3; MAPKKKK3; RAB8IPL1;

**Catalog #:** A4224

**Antibody Type:**

Monoclonal Antibody

**Species:** 23

**Gene ID:** 8491

**Isotype:** IgG

**Swiss Prot:** Q8IVH8

**Purity:** Affinity purification

For research use only.

**Background:**

This gene encodes a member of the mitogen-activated protein kinase kinase kinase family.

The encoded protein activates key effectors in cell signalling, among them c-Jun. Alternatively spliced transcripts encoding multiple isoforms have been observed for this gene.

*To place an order, please [Click HERE](#).*