

ASMTL

Reactivity:Human

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:20 - 1:200 IF 1:20 - 1:100

Calculated MW:69kDa

Observed MW:Refer to figures

Immunogen:

A synthetic Peptide of human ASMTL

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

f

Synonym:

ASTML; ASMTLX; ASMTLY;

Catalog #:A4237

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:8623

Isotype:IgG

Swiss Prot:O95671

Purity:Affinity purification

For research use only.

Background:

The protein encoded by this gene has an N-terminus that is similar to the multicopy associated filamentation (maf) protein of *Bacillus subtilis* and to orfE of *Escherichia coli*, while the C-terminus is similar to N-acetylserotonin O-methyltransferase. This gene is located in the pseudoautosomal region 1 (PAR1) of X and Y chromosomes. Three transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).