

SSNA1

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:20 - 1:200

Calculated MW:14kDa

Observed MW:Refer to figures

Immunogen:

A synthetic Peptide of human SSNA1

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

d

Synonym:

N14; NA14; NA-14;

To place an order, please [Click HERE](#).

Catalog #:A4239

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:8636

Isotype:IgG

Swiss Prot:O43805

Purity:Affinity purification

For research use only.