

EIF3G

Reactivity:Human Mouse Rat

Tested applications:WB

Recommended Dilution:WB 1:200 - 1:2000

Calculated MW:36kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human EIF3G

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

j

Synonym:

EIF3S4; EIF3-P42; eIF3-p44; eIF3-delta;

To place an order, please [Click HERE](#).

Catalog #:A4240

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:8666

Isotype:IgG

Swiss Prot:O75821

Purity:Affinity purification

For research use only.