

SNX3

Reactivity:Human Mouse Rat

Tested applications:WB IF

Recommended Dilution:WB 1:200 - 1:2000 IF 1:10 - 1:100

Calculated MW:19kDa

Observed MW:Refer to figures

Immunogen:

A synthetic Peptide of human SNX3

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

SDP3; Grd19; MCOPS8;

Catalog #:A4246

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:8724

Isotype:IgG

Swiss Prot:O60493

Purity:Affinity purification

For research use only.

Background:

This gene encodes a member of the sorting nexin family. Members of this family contain a phox (PX) domain, which is a phosphoinositide binding domain, and are involved in intracellular trafficking. This protein does not contain a coiled coil region, like most family members. This protein interacts with phosphatidylinositol-3-phosphate, and is involved in protein trafficking. A pseudogene of this gene is present on the sex chromosomes. Alternative splicing results in multiple transcript variants encoding distinct isoforms.

To place an order, please [Click HERE](#).