

URI1

Reactivity:Human

Tested applications:WB IHC IP

Recommended Dilution:WB 1:200 - 1:2000 IHC 1:20 - 1:200 IP 1:20 - 1:50

Calculated MW:60kDa

Observed MW:Refer to figures

Immunogen:

A synthetic Peptide of human URI1

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

RMP; URI; NNX3; C19orf2; PPP1R19;

Catalog #:A4247

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:8725

Isotype:IgG

Swiss Prot:O94763

Purity:Affinity purification

For research use only.

Background:

This gene encodes member of the prefoldin family of molecular chaperones. The encoded protein functions as a scaffolding protein and plays roles in ubiquitination and transcription, in part through interactions with the RNA polymerase II subunit RPB5. This gene may play a role in multiple malignancies including ovarian cancer and hepatocellular carcinoma. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene, and a pseudogene of this gene is located on the long arm of chromosome 22.

To place an order, please [Click HERE](#).