

PSTPIP1

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:200 - 1:1000 IHC 1:20 - 1:200

Calculated MW:48kDa

Observed MW:Refer to figures

Immunogen:

A synthetic Peptide of human PSTPIP1

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

H-PIP; PAPAS; CD2BP1; PSTPIP; CD2BP1L; CD2BP1S;

Catalog #:A4266

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:9051

Isotype:IgG

Swiss Prot:O43586

Purity:Affinity purification

For research use only.

Background:

The protein encoded by this gene binds to the cytoplasmic tail of CD2, an effector of T cell activation and adhesion, negatively affecting CD2-triggered T cell activation. The encoded protein appears to be a scaffold protein and a regulator of the actin cytoskeleton. It has also been shown to bind ABL1, PTPN18, WAS, CD2AP, and PTPN12. Mutations in this gene are a cause of PAPA syndrome.

To place an order, please [Click HERE](#).