

EIF1AY

Reactivity:Human Mouse Rat

Tested applications:WB

Recommended Dilution:WB 1:200 - 1:2000

Calculated MW:16kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human EIF1AY

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

eIF-4C;

Catalog #:A4270

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:9086

Isotype:IgG

Swiss Prot:O14602

Purity:Affinity purification

For research use only.

Background:

This gene is located on the non-recombining region of the Y chromosome. It encodes a protein related to eukaryotic translation initiation factor 1A (EIF1A), which may function in stabilizing the binding of the initiator Met-tRNA to 40S ribosomal subunits. Alternative splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).