

LONP1

Reactivity:Human Mouse Rat

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:106kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human LONP1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

h

Synonym:

LON; LONP; PIM1; hLON; LonHS; CODASS; PRSS15;

Catalog #:A4293

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:9361

Isotype:IgG

Swiss Prot:P36776

Purity:Affinity purification

For research use only.

Background:

This gene encodes a mitochondrial matrix protein that belongs to the Lon family of ATP-dependent proteases. This protein mediates the selective degradation of misfolded, unassembled or oxidatively damaged polypeptides in the mitochondrial matrix. It may also have a chaperone function in the assembly of inner membrane protein complexes, and participate in the regulation of mitochondrial gene expression and maintenance of the integrity of the mitochondrial genome. Decreased expression of this gene has been noted in a patient with hereditary spastic paraplegia (PMID:18378094). Alternatively spliced transcript variants have been found for this gene.

To place an order, please [Click HERE](#).