

PLAA

Reactivity:Human Mouse

Tested applications:WB IP

Recommended Dilution:WB 1:200 - 1:2000 IP 1:20 - 1:100

Calculated MW:87kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human PLAA

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

DOA1; PLAP; PLA2P;

To place an order, please [Click HERE](#).

Catalog #:A4294

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:9373

Isotype:IgG

Swiss Prot:Q9Y263

Purity:Affinity purification

For research use only.