

GRHPR

Reactivity:Human Mouse

Tested applications:WB IHC

Recommended Dilution:WB 1:200 - 1:2000 IHC 1:20 - 1:200

Calculated MW:36kDa

Observed MW:Refer to figures

Immunogen:

A synthetic Peptide of human GRHPR

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

PH2; GLXR; GLYD;

Catalog #:A4295

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:9380

Isotype:IgG

Swiss Prot:Q9UBQ7

Purity:Affinity purification

For research use only.

Background:

This gene encodes an enzyme with hydroxypyruvate reductase, glyoxylate reductase, and D-glycerate dehydrogenase enzymatic activities. The enzyme has widespread tissue expression and has a role in metabolism. Type II hyperoxaluria is caused by mutations in this gene.

To place an order, please [Click HERE](#).