

RPH3AL

Reactivity:Human

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW:34kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human RPH3AL

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

NOC2;

Catalog #:A4311

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:9501

Isotype:IgG

Swiss Prot:Q9UNE2

Purity:Affinity purification

For research use only.

Background:

The protein encoded by this gene plays a direct regulatory role in calcium-ion-dependent exocytosis in both endocrine and exocrine cells and plays a key role in insulin secretion by pancreatic cells. This gene is likely a tumor suppressor. Alternative splicing results in multiple transcript variants encoding distinct isoforms.

To place an order, please [Click HERE](#).