

PMPCB

Reactivity:Human Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:54kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human PMPCB

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

h

Synonym:

MPPB; P-52; MPP11; MPPP52; Beta-MPP;

Catalog #:A4312

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:9512

Isotype:IgG

Swiss Prot:O75439

Purity:Affinity purification

For research use only.

Background:

This gene is a member of the peptidase M16 family and encodes a protein with a zinc-binding motif. This protein is located in the mitochondrial matrix and catalyzes the cleavage of the leader peptides of precursor proteins newly imported into the mitochondria, though it only functions as part of a heterodimeric complex.

To place an order, please [Click HERE](#).