

## EEF1E1

**Reactivity:**Human Mouse

**Tested applications:**WB FC

**Recommended Dilution:**WB 1:200 - 1:2000 FC 1:20 - 1:50

**Calculated MW:**20kDa

**Observed MW:**Refer to figures

**Immunogen:**

A synthetic peptide of human EEF1E1

**Storage Buffer:**

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Concentration:**

b

**Synonym:**

P18; AIMP3;

**Catalog #:**A4314

**Antibody Type:**

Monoclonal Antibody

**Species:**Mouse

**Gene ID:**9521

**Isotype:**IgG

**Swiss Prot:**O43324

**Purity:**Affinity purification

For research use only.

**Background:**

This gene encodes a multifunctional protein that localizes to both the cytoplasm and nucleus. In the cytoplasm, the encoded protein is an auxiliary component of the macromolecular aminoacyl-tRNA synthase complex. However, its mouse homolog has been shown to translocate to the nucleus in response to DNA damage, and it plays a positive role in ATM/ATR-mediated p53 activation. Alternative splicing results in multiple transcript variants. Read-through transcription also exists between this gene and the neighboring downstream MUTED (muted homolog) gene. An EEF1E1-related pseudogene has been identified on chromosome 2.

**To place an order, please [Click HERE](#).**