

SEC22B

Reactivity:Human

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:24kDa

Observed MW:Refer to Figures

Immunogen:

A synthetic peptide of human SEC22B

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ERS-24; SEC22L1

Catalog #:A4318

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:9554

Isotype:IgG

Swiss Prot:O75396

Purity:Affinity purification

For research use only.

Background:

The protein encoded by this gene is a member of the SEC22 family of vesicle trafficking proteins. It seems to complex with SNARE and it is thought to play a role in the ER-Golgi protein trafficking. This protein has strong similarity to Mus musculus and Cricetulus griseus proteins.

To place an order, please [Click HERE](#).