

SEC22B

Reactivity: Human

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 24kDa

Observed MW: Refer to Figures

Immunogen:

A synthetic peptide of human SEC22B

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ERS-24; SEC22L1

Catalog #: A4318

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 9554

Isotype: IgG

Swiss Prot: O75396

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a member of the SEC22 family of vesicle trafficking proteins. It seems to complex with SNARE and it is thought to play a role in the ER-Golgi protein trafficking. This protein has strong similarity to Mus musculus and Cricetulus griseus proteins.

To place an order, please [Click HERE](#).