

C14orf2

Reactivity: Human

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:20 - 1:200

Calculated MW: 7kDa

Observed MW: Refer to figures

Immunogen:

A synthetic Peptide of human C14orf2

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

MP68; PLPM;

To place an order, please [Click HERE](#).

Catalog #: A4319

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 9556

Isotype: IgG

Swiss Prot: P56378

Purity: Affinity purification

For research use only.