

## GOSR2

**Reactivity:** Human

**Tested applications:** WB

**Recommended Dilution:** WB 1:500 - 1:2000

**Calculated MW:** 25kDa

**Observed MW:** Refer to figures

**Immunogen:**

Recombinant protein of human GOSR2

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

Bos1; EPM6; GS27;

**Catalog #:** A4321

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 9570

**Isotype:** IgG

**Swiss Prot:** O14653

**Purity:** Affinity purification

For research use only.

**Background:**

This gene encodes a trafficking membrane protein which transports proteins among the medial- and trans-Golgi compartments. Due to its chromosomal location and trafficking function, this gene may be involved in familial essential hypertension. Three transcript variants encoding three different isoforms have been found for this gene.

*To place an order, please [Click HERE](#).*