

GOSR2

Reactivity:Human

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:25kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human GOSR2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

Bos1; EPM6; GS27;

Catalog #:A4321

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:9570

Isotype:IgG

Swiss Prot:O14653

Purity:Affinity purification

For research use only.

Background:

This gene encodes a trafficking membrane protein which transports proteins among the medial- and trans-Golgi compartments. Due to its chromosomal location and trafficking function, this gene may be involved in familial essential hypertension. Three transcript variants encoding three different isoforms have been found for this gene.

To place an order, please [Click HERE](#).