

## PREPL

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB

**Recommended Dilution:**WB 1:500 - 1:2000

**Calculated MW:**84kDa

**Observed MW:**Refer to figures

**Immunogen:**

A synthetic Peptide of human PREPL

**Storage Buffer:**

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Catalog #:**A4322

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**9581

**Isotype:**IgG

**Swiss Prot:**Q4J6C6

**Purity:**Affinity purification

For research use only.

**Background:**

The protein encoded by this gene belongs to the prolyl oligopeptidase subfamily of serine peptidases. Mutations in this gene have been associated with hypotonia-cystinuria syndrome, also known as the 2p21 deletion syndrome. Several alternatively spliced transcript variants encoding either the same or different isoforms have been described for this gene.

**To place an order, please [Click HERE](#).**