

PDIA4

Reactivity:Human

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:73kDa

Observed MW:Refer to Figures

Immunogen:

A synthetic peptide of human PDIA4

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

j

Synonym:

ERP70; ERP72; ERp-72;

To place an order, please [Click HERE](#).

Catalog #:A4326

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:9601

Isotype:IgG

Swiss Prot:P13667

Purity:Affinity purification

For research use only.