

PIEZO1

Reactivity:Human

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:287kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human PIEZO1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

DHS; Mib; FAM38A;

Catalog #:A4340

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:9780

Isotype:IgG

Swiss Prot:Q92508

Purity:Affinity purification

For research use only.

Background:

The protein encoded by this gene is a mechanically-activated ion channel that links mechanical forces to biological signals. The encoded protein contains 36 transmembrane domains and functions as a homotetramer. Defects in this gene have been associated with dehydrated hereditary stomatocytosis.

To place an order, please [Click HERE](#).