

## SRGAP3

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB IHC IP

**Recommended Dilution:**WB 1:500 - 1:2000 IHC 1:20 - 1:200 IP 1:20 - 1:50

**Calculated MW:**125kDa

**Observed MW:**Refer to figures

**Immunogen:**

A synthetic Peptide of human SRGAP3

**Storage Buffer:**

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

WRP; MEGAP; SRGAP2; ARHGAP14;

**To place an order, please [Click HERE](#).**

**Catalog #:**A4354

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**9901

**Isotype:**IgG

**Swiss Prot:**O43295

**Purity:**Affinity purification

For research use only.