

## GNPDA1

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB IHC

**Recommended Dilution:**WB 1:500 - 1:2000 IHC 1:20 - 1:200

**Calculated MW:**33kDa

**Observed MW:**Refer to figures

**Immunogen:**

A synthetic Peptide of human GNPDA1

**Storage Buffer:**

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

GPI; HLN; GNP1; GNPI; GNPDA;

**Catalog #:**A4358

**Antibody Type:**

Monoclonal Antibody

**Species:**23

**Gene ID:**10007

**Isotype:**IgG

**Swiss Prot:**P46926

**Purity:**Affinity purification

For research use only.

**Background:**

Glucosamine-6-phosphate deaminase (EC 3.5.99.6) is an allosteric enzyme that catalyzes the reversible conversion of D-glucosamine-6-phosphate into D-fructose-6-phosphate and ammonium (Arreola et al., 2003)

*To place an order, please [Click HERE](#).*