

## TNIP2

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB IHC IF

**Recommended Dilution:**WB 1:500 - 1:2000 IHC 1:20 - 1:200 IF 1:50 - 1:100

**Calculated MW:**49kDa

**Observed MW:**Refer to figures

**Immunogen:**

Recombinant protein of human TNIP2

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

KLIP; ABIN2; FLIP1;

**Catalog #:**A4962

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**79155

**Isotype:**IgG

**Swiss Prot:**Q8NFZ5

**Purity:**Affinity purification

For research use only.

**Background:**

This gene encodes a protein which acts as an inhibitor of NFkappaB activation. The encoded protein is also involved in MAP/ERK signaling pathway in specific cell types. It may be involved in apoptosis of endothelial cells. Alternative splicing results in multiple transcript variants. A pseudogene related to this gene is located on the X chromosome.

*To place an order, please [Click HERE](#).*