

CHMP6

Reactivity:Human Mouse Rat

Tested applications:WB IP

Recommended Dilution:WB 1:200 - 1:2000 IP 1:20 - 1:50

Calculated MW:23kDa

Observed MW:Refer to figures

Immunogen:

A synthetic Peptide of human CHMP6

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

VPS20;

Catalog #:A4975

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:79643

Isotype:IgG

Swiss Prot:Q96FZ7

Purity:Affinity purification

For research use only.

Background:

This gene encodes a member of the chromatin-modifying protein/charged multivesicular body protein family. Proteins in this family are part of the ESCRT-III (endosomal sorting complex required for transport III) which degrades surface receptors, and in biosynthesis of endosomes.

To place an order, please [Click HERE](#).