

GTDC1

Reactivity:Human Mouse Rat

Tested applications:WB IHC IF

Recommended Dilution:WB 1:200 - 1:2000 IHC 1:20 - 1:200 IF 1:10 - 1:100

Calculated MW:53kDa

Observed MW:Refer to figures

Immunogen:

A synthetic Peptide of human GTDC1

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

mat-Xa; Hmat-Xa;

To place an order, please [Click HERE](#).

Catalog #:A4980

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:79712

Isotype:IgG

Swiss Prot:Q4AE62

Purity:Affinity purification

For research use only.