

NHEJ1

Reactivity:Human Mouse

Tested applications:WB

Recommended Dilution:WB 1:200 - 1:2000

Calculated MW:33kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human NHEJ1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

e

Synonym:

XLf;

Catalog #:A4985

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:79840

Isotype:IgG

Swiss Prot:Q9H9Q4

Purity:Affinity purification

For research use only.

Background:

Double-strand breaks in DNA result from genotoxic stresses and are among the most damaging of DNA lesions. This gene encodes a DNA repair factor essential for the nonhomologous end-joining pathway, which preferentially mediates repair of double-stranded breaks. Mutations in this gene cause different kinds of severe combined immunodeficiency disorders.

To place an order, please [Click HERE](#).