

LPCAT1

Reactivity:Human Mouse

Tested applications:WB

Recommended Dilution:WB 1:1000 - 1:2000

Calculated MW:59kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human LPCAT1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

AYTL2; lpcat; PFAAP3;

Background:

Lysophosphatidylcholine (LPC) acyltransferase (LPCAT; EC 2.3.1.23) catalyzes the conversion of LPC to phosphatidylcholine (PC) in the remodeling pathway of PC biosynthesis (Nakanishi et al., 2006

To place an order, please [Click HERE](#).

Catalog #:A4987

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:79888

Isotype:IgG

Swiss Prot:Q8NF37

Purity:Affinity purification

For research use only.