

DNAJB14

Reactivity:Human Mouse Rat

Tested applications:WB

Recommended Dilution:WB 1:1000 - 1:4000

Calculated MW:43kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human DNAJB14

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

f

Synonym:

EGNR9427; PRO34683;

To place an order, please [Click HERE](#).

Catalog #:A4990

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:79982

Isotype:IgG

Swiss Prot:Q8TBM8

Purity:Affinity purification

For research use only.