

PTGES2

Reactivity:Human Mouse Rat

Tested applications:WB IHC IF IP

Recommended Dilution:WB 1:200 - 1:2000 IHC 1:20 - 1:200 IF 1:10 - 1:100 IP 1:20 - 1:50

Calculated MW:42kDa

Observed MW:Refer to figures

Immunogen:

A synthetic Peptide of human PTGES2

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

GBF1; GBF-1; PGES2; C9orf15; mPGES-2;

Catalog #:A4993

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:80142

Isotype:IgG

Swiss Prot:Q9H7Z7

Purity:Affinity purification

For research use only.

Background:

The protein encoded by this gene is a membrane-associated prostaglandin E synthase, which catalyzes the conversion of prostaglandin H2 to prostaglandin E2. This protein also has been shown to activate the transcription regulated by a gamma-interferon-activated transcription element (GATE). Multiple transcript variants have been found for this gene.

To place an order, please [Click HERE](#).