

EDC3

Reactivity: Human Mouse Rat

Tested applications: WB IHC IP

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:20 - 1:200 IP 1:20 - 1:50

Calculated MW: 56kDa

Observed MW: Refer to figures

Immunogen:

A synthetic Peptide of human EDC3

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

YJDC; LSM16; YJEFN2;

Catalog #: A4994

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 80153

Isotype: IgG

Swiss Prot: Q96F86

Purity: Affinity purification

For research use only.

Background:

EDC3 is associated with an mRNA-decapping complex required for removal of the 5-prime cap from mRNA prior to its degradation from the 5-prime end (Fenger-Gron et al., 2005)

To place an order, please [Click HERE](#).