

GRPEL1

Reactivity:Human Mouse Rat

Tested applications:WB IHC IF IP

Recommended Dilution:WB 1:200 - 1:2000 IHC 1:20 - 1:200 IF 1:10 - 1:100 IP 1:20 - 1:50

Calculated MW:24kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human GRPEL1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

h

Synonym:

GrpE; HMGE;

To place an order, please [Click HERE](#).

Catalog #:A4999

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:80273

Isotype:IgG

Swiss Prot:Q9HAV7

Purity:Affinity purification

For research use only.