

ILKAP

Reactivity:Human Mouse Rat

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:43kDa

Observed MW:Refer to figures

Immunogen:

A synthetic Peptide of human ILKAP

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ILKAP2; ILKAP3; PP2C-DELTA;

Catalog #:A5010

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:80895

Isotype:IgG

Swiss Prot:Q9H0C8

Purity:Affinity purification

For research use only.

Background:

The protein encoded by this gene is a protein serine/threonine phosphatase of the PP2C family.

This protein can interact with integrin-linked kinase (ILK/ILK1), a regulator of integrin mediated signaling, and regulate the kinase activity of ILK. Through the interaction with ILK, this protein may selectively affect the signaling process of ILK-mediated glycogen synthase kinase 3 beta (GSK3beta), and thus participate in Wnt signaling pathway.

To place an order, please [Click HERE](#).