

ZBP1

Reactivity:Human

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:20 - 1:200

Calculated MW:46kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human ZBP1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

DAI; DLM1; DLM-1; C20orf183;

Background:

This gene encodes a Z-DNA binding protein. The encoded protein plays a role in the innate immune response by binding to foreign DNA and inducing type-I interferon production.

Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene.

To place an order, please [Click HERE](#).

Catalog #:A5011

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:81030

Isotype:IgG

Swiss Prot:Q9H171

Purity:Affinity purification

For research use only.