

## DIAPH3

---

**Reactivity:**Human

**Tested applications:**WB IHC IF

**Recommended Dilution:**WB 1:500 - 1:2000 IHC 1:20 - 1:200 IF 1:20 - 1:100

**Calculated MW:**137kDa

**Observed MW:**Refer to figures

**Immunogen:**

A synthetic Peptide of human DIAPH3

**Storage Buffer:**

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Concentration:**

r

**Synonym:**

AN; DIA2; DRF3; AUNA1; NSDAN; diap3; mDia2;

**Catalog #:**A5019

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**81624

**Isotype:**IgG

**Swiss Prot:**Q9NSV4

**Purity:**Affinity purification

For research use only.

**Background:**

This gene encodes a member of the diaphanous subfamily of the formin family. Members of this family are involved in actin remodeling and regulate cell movement and adhesion. Mutations in this gene are associated with autosomal dominant auditory neuropathy 1. Multiple transcript variants encoding different isoforms have been found for this gene.

*To place an order, please [Click HERE](#).*