

BCO2

Reactivity:Human Mouse Rat

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:66kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human BCO2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

BCDO2; B-DIOX-II;

Catalog #:A5034

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:83875

Isotype:IgG

Swiss Prot:Q9BYV7

Purity:Affinity purification

For research use only.

Background:

This gene encodes an enzyme which oxidizes carotenoids such as beta-carotene during the biosynthesis of vitamin A. Multiple transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).