

POLDIP3

Reactivity:Human Mouse Rat

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:46kDa

Observed MW:Refer to figures

Immunogen:

A synthetic Peptide of human POLDIP3

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

j

Synonym:

SKAR; PDIP46;

Catalog #:A5048

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:84271

Isotype:IgG

Swiss Prot:Q9BY77

Purity:Affinity purification

For research use only.

Background:

This gene encodes an RRM (RNA recognition motif)-containing protein that participates in the regulation of translation by recruiting ribosomal protein S6 kinase beta-1 to mRNAs. Alternative splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).