

SARNP

Reactivity:Human Mouse Rat

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:24kDa

Observed MW:Refer to figures

Immunogen:

A synthetic Peptide of human SARNP

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

HCC1; THO1; CIP29; HSPC316;

Catalog #:A5054

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:84324

Isotype:IgG

Swiss Prot:P82979

Purity:Affinity purification

For research use only.

Background:

This gene encodes a protein that is upregulated in response to various cytokines. The encoded protein may play a role in cell cycle progression. A translocation between this gene and the myeloid/lymphoid leukemia gene, resulting in expression of a chimeric protein, has been associated with acute myelomonocytic leukemia. Pseudogenes exist on chromosomes 7 and 8. Alternatively spliced transcript variants have been described.

To place an order, please [Click HERE](#).