

ACCS

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:20 - 1:200

Calculated MW:57kDa

Observed MW:Refer to figures

Immunogen:

A synthetic Peptide of human ACCS

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

f

Synonym:

ACS; PHACS;

To place an order, please [Click HERE](#).

Catalog #:A5062

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:84680

Isotype:IgG

Swiss Prot:Q96QU6

Purity:Affinity purification

For research use only.