

CBR4

Reactivity:Human Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:200 - 1:2000 IHC 1:20 - 1:200

Calculated MW:25kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human CBR4

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

SDR45C1;

To place an order, please [Click HERE](#).

Catalog #:A5069

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:84869

Isotype:IgG

Swiss Prot:Q8N4T8

Purity:Affinity purification

For research use only.