

CALML4

Reactivity:Human Mouse Rat

Tested applications:WB IF IP

Recommended Dilution:WB 1:200 - 1:2000 IF 1:10 - 1:100 IP 1:20 - 1:50

Calculated MW:22kDa

Observed MW:Refer to figures

Immunogen:

A synthetic Peptide of human CALML4

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

NY-BR-20;

To place an order, please [Click HERE](#).

Catalog #:A5086

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:91860

Isotype:IgG

Swiss Prot:Q96GE6

Purity:Affinity purification

For research use only.