

MOB1B

Reactivity:Human Mouse Rat

Tested applications:WB IHC IP

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:20 - 1:200 IP 1:20 - 1:50

Calculated MW:25kDa

Observed MW:Refer to figures

Immunogen:

A synthetic Peptide of human MOB1B

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

MATS2; MOB4A; MOBKL1A;

Catalog #:A5091

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:92597

Isotype:IgG

Swiss Prot:Q7L9L4

Purity:Affinity purification

For research use only.

Background:

The protein encoded by this gene is similar to the yeast Mob1 protein. Yeast Mob1 binds Mps1p, a protein kinase essential for spindle pole body duplication and mitotic checkpoint regulation. Three transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).