

ATPIF1

Reactivity:Human Mouse

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW:12kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human ATPIF1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

a

Synonym:

IP; ATPi; ATPiP;

Catalog #:A5099

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:93974

Isotype:IgG

Swiss Prot:Q9UII2

Purity:Affinity purification

For research use only.

Background:

This gene encodes a mitochondrial ATPase inhibitor. Alternative splicing occurs at this locus and three transcript variants encoding distinct isoforms have been identified.

To place an order, please [Click HERE](#).