

## MYL12B

---

**Reactivity:**Mouse Rat

**Tested applications:**WB

**Recommended Dilution:**WB 1:500 - 1:2000

**Calculated MW:**20kDa

**Observed MW:**Refer to figures

**Immunogen:**

Recombinant protein of human MYL12B

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

MLC-B; MRLC2;

**Background:**

The activity of nonmuscle myosin II (see MYH9; MIM 160775) is regulated by phosphorylation of a regulatory light chain, such as MRLC2. This phosphorylation results in higher MgATPase activity and the assembly of myosin II filaments (Iwasaki et al., 2001)

*To place an order, please [Click HERE](#).*

**Catalog #:**A5104

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**103910

**Isotype:**IgG

**Swiss Prot:**O14950

**Purity:**Affinity purification

For research use only.