

PIK3IP1

Reactivity:Human Mouse Rat

Tested applications:WB IP

Recommended Dilution:WB 1:500 - 1:2000 IP 1:20 - 1:50

Calculated MW:28kDa

Observed MW:Refer to figures

Immunogen:

A synthetic Peptide of human PIK3IP1

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

HGFL; hHGFL(S);

To place an order, please [Click HERE](#).

Catalog #:A5110

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:113791

Isotype:IgG

Swiss Prot:Q96FE7

Purity:Affinity purification

For research use only.