

## KCTD12

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB IF

**Recommended Dilution:**WB 1:500 - 1:2000 IF 1:10 - 1:100

**Calculated MW:**36kDa

**Observed MW:**Refer to figures

**Immunogen:**

A synthetic Peptide of human KCTD12

**Storage Buffer:**

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Concentration:**

g

**Synonym:**

PFET1; PFETIN; C13orf2;

*To place an order, please [Click HERE](#).*

**Catalog #:**A5116

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**115207

**Isotype:**IgG

**Swiss Prot:**Q96CX2

**Purity:**Affinity purification

For research use only.